

RIOK3 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP01940-1**

Capture Antibody Information

Catalog Number: 85449-3-PBS	Clone ID: 242996F11	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: RIO kinase 3 (yeast)
Isotype: IgG	GenBank: BC039729	Gene ID: 8780
Purification Method: Protein A purification	Immunogen Catalog Number: Ag4497	

Detection Antibody Information

Catalog Number: 85449-1-PBS	Clone ID: 242996C6	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: RIO kinase 3 (yeast)
Isotype: IgG	GenBank: BC039729	Gene ID: 8780
Purification Method: Protein A purification	Immunogen Catalog Number: Ag4497	

Applications

Tested Applications: Cytometric bead array	Range: 0.781-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

Product Information

MP01940-1 targets RIOK3 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: RIOK3 Recombinant antibody, PBS Only (Capture) 85449-3-PBS (242996F11). 100 µg. Concentration 1 mg/mL.

Detection antibody: RIOK3 Recombinant antibody, PBS Only (Detector) 85449-1-PBS (242996C6). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP01940-1, R1OK3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85449-3-PBS. Detection antibody: 85449-1-PBS. Standard: Ag4497. Range: 0.781-100 ng/mL.