

Mouse IL-28A Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP01988-2**

Capture Antibody Information

Catalog Number:

98388-1-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

242247E2

Reactivity:

mouse

GenBank:

NM_001024673.2

Conjugate:

Unconjugated

Full name:

interleukin 28A

Gene ID:

330496

Detection Antibody Information

Catalog Number:

85590-3-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

242247B2

Reactivity:

mouse

GenBank:

NM_001024673.2

Conjugate:

Unconjugated

Full name:

interleukin 28A

Gene ID:

330496

Applications

Tested Applications:

Cytometric bead array

Range:

0.313-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:

It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP01988-2 targets IL-28A in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Anti-Mouse IL28a Rabbit Recombinant Antibody, PBS Only (Capture) 98388-1-PBS (242247E2). 100 µg. Concentration 1 mg/mL.

Detection antibody: Mouse IL28a Recombinant antibody, PBS Only (Capture/Detector) 85590-3-PBS (242247B2). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:


Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:

PBS only

Selected Validation Data



Cytometric bead array standard curve of MP01988-2, MOUSE IL-28A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 98388-1-PBS. Detection antibody: 85590-3-PBS. Standard: Eg3169. Range: 0.313-20 ng/mL