

SF3B4 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP02012-1

Capture Antibody Information

Catalog Number:
85663-3-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
243098E9

Reactivity:
human

GenBank:
BC004273

Immunogen Catalog Number:
Ag0737

Conjugate:
Unconjugated

Full name:
splicing factor 3b, subunit 4, 49kDa

Gene ID:
10262

Detection Antibody Information

Catalog Number:
85663-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
243098B6

Reactivity:
human

GenBank:
BC004273

Immunogen Catalog Number:
Ag0737

Conjugate:
Unconjugated

Full name:
splicing factor 3b, subunit 4, 49kDa

Gene ID:
10262

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP02012-1 targets SF3B4 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: SF3B4 Recombinant antibody, PBS Only (Capture/Detector) 85663-3-PBS (243098E9). 100 µg. Concentration 1 mg/mL.

Detection antibody: SF3B4 Recombinant antibody, PBS Only (Detector) 85663-2-PBS (243098B6). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP02012-1, SF3B4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85663-3-PBS. Detection antibody: 85663-2-PBS. Standard: Ag0737. Range: 0.781-100 ng/mL.