

RNF20 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP02014-1

Capture Antibody Information

Catalog Number:	Clone ID:	Conjugate:
85630-6-PBS	243078G12	Unconjugated
Host:	Reactivity:	Full name:
Rabbit	human, mouse	ring finger protein 20
Isotype:	GenBank:	Gene ID:
IgG	BC110585	56254
Purification Method:	Immunogen Catalog Number:	
Protein A purification	Ag16066	

Detection Antibody Information

Catalog Number:	Clone ID:	Conjugate:
85630-5-PBS	243078F6	Unconjugated
Host:	Reactivity:	Full name:
Rabbit	human	ring finger protein 20
Isotype:	GenBank:	Gene ID:
IgG	BC110585	56254
Purification Method:	Immunogen Catalog Number:	
Protein A purification	Ag16066	

Applications

Tested Applications:	Range:	Recommended Dilutions:
Cytometric bead array	0.781-100 ng/mL (Cytometric Bead Array)	It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP02014-1 targets RNF20 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: RNF20 Recombinant antibody, PBS Only (Capture) 85630-6-PBS (243078G12). 100 μ g. Concentration 1 mg/mL.

Detection antibody: RNF20 Recombinant antibody, PBS Only (Capture/Detector) 85630-5-PBS (243078F6). 100 μ g. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP02014-1, RNF20 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85630-6-PBS. Detection antibody: 85630-5-PBS. Standard: Ag16066. Range: 0.781-100 ng/mL.