

KTN1 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP02037-2**

Capture Antibody Information

Catalog Number: 85681-3-PBS	Clone ID: 243067D10	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: kinectin 1 (kinesin receptor)
Isotype: IgG	GenBank: BC117132	Gene ID: 3895
Purification Method: Protein A purification	Immunogen Catalog Number: Ag13854	

Detection Antibody Information

Catalog Number: 85681-1-PBS	Clone ID: 243067F10	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: kinectin 1 (kinesin receptor)
Isotype: IgG	GenBank: BC117132	Gene ID: 3895
Purification Method: Protein A purification	Immunogen Catalog Number: Ag13854	

Applications

Tested Applications: Cytometric bead array	Range: 3.125-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

Product Information

MP02037-2 targets KTN1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: KTN1 Recombinant antibody, PBS Only (Capture) 85681-3-PBS (243067D10). 100 µg. Concentration 1 mg/mL.

Detection antibody: KTN1 Recombinant antibody, PBS Only (Detector) 85681-1-PBS (243067F10). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP02037-2, KTN1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85681-3-PBS. Detection antibody: 85681-1-PBS. Standard: Ag13854. Range: 3.125-100 ng/mL.