

STIM2 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP02057-1**

Capture Antibody Information

Catalog Number: 85727-5-PBS	Clone ID: 250016D4	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: stromal interaction molecule 2
Isotype: IgG	GenBank: BC146661	Gene ID: 57620
Purification Method: Protein A purification	Immunogen Catalog Number: Ag15530	

Detection Antibody Information

Catalog Number: 85727-4-PBS	Clone ID: 250016F10	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: stromal interaction molecule 2
Isotype: IgG	GenBank: BC146661	Gene ID: 57620
Purification Method: Protein A purification	Immunogen Catalog Number: Ag15530	

Applications

Tested Applications: Cytometric bead array	Range: 0.781-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

Product Information

MP02057-1 targets STIM2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: STIM2 Recombinant antibody, PBS Only (Capture) 85727-5-PBS (250016D4). 100 µg. Concentration 1 mg/mL.

Detection antibody: STIM2 Recombinant antibody, PBS Only (Detector) 85727-4-PBS (250016F10). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP02057-1, STIM2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85727-5-PBS. Detection antibody: 85727-4-PBS. Standard: Ag15530. Range: 0.781-100 ng/mL.