

# Centrin 1 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP02059-2**

## Capture Antibody Information

Catalog Number: 85637-2-PBS	Clone ID: 243013C4	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: centrin, EF-hand protein, 1
Isotype: IgG	GenBank: BC029515	Gene ID: 1068
Purification Method: Protein A purification	Immunogen Catalog Number: Ag3529	

## Detection Antibody Information

Catalog Number: 85637-1-PBS	Clone ID: 243013C11	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: centrin, EF-hand protein, 1
Isotype: IgG	GenBank: BC029515	Gene ID: 1068
Purification Method: Protein A purification	Immunogen Catalog Number: Ag3529	

## Applications

Tested Applications: Cytometric bead array	Range: 0.781-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
-----------------------------------------------	---------------------------------------------------	------------------------------------------------------------------------------------------------------------------------------------

## Product Information

MP02059-2 targets Centrin 1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CETN1 Recombinant antibody, PBS Only (Capture) 85637-2-PBS (243013C4). 100 µg. Concentration 1 mg/mL.

Detection antibody: CETN1 Recombinant antibody, PBS Only (Detector) 85637-1-PBS (243013C11). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
Storage buffer:  
PBS only

# Selected Validation Data



Cytometric bead array standard curve of MP02059-2, Centrin 1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85637-2-PBS. Detection antibody: 85637-1-PBS. Standard: Ag3529. Range: 0.781-100 ng/mL.