

FARP1 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP02060-1**

Capture Antibody Information

Catalog Number:
85638-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
243132A3
Reactivity:
human
GenBank:
BC071592
Immunogen Catalog Number:
Ag21225

Conjugate:
Unconjugated
Full name:
FERM, RhoGEF (ARHGEF) and pleckstrin domain protein 1 (chondrocyte-derived)
Gene ID:
10160

Detection Antibody Information

Catalog Number:
85638-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
243132B6
Reactivity:
human
GenBank:
BC071592
Immunogen Catalog Number:
Ag21225

Conjugate:
Unconjugated
Full name:
FERM, RhoGEF (ARHGEF) and pleckstrin domain protein 1 (chondrocyte-derived)
Gene ID:
10160

Applications

Tested Applications:
Cytometric bead array

Range:
0.469-60 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP02060-1 targets FARP1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: FARP1 Recombinant antibody, PBS Only (Capture) 85638-3-PBS (243132A3). 100 µg. Concentration 1 mg/mL.

Detection antibody: FARP1 Recombinant antibody, PBS Only (Detector) 85638-2-PBS (243132B6). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP02060-1, FARP1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85638-3-PBS. Detection antibody: 85638-2-PBS. Standard: Ag21225. Range: 0.469-60 ng/mL