

Neurochondrin Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP02072-2**

Capture Antibody Information

Catalog Number:
85730-4-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
243146D9
Reactivity:
human
GenBank:
BC024592
Immunogen Catalog Number:
Ag3859

Conjugate:
Unconjugated
Full name:
neurochondrin
Gene ID:
23154

Detection Antibody Information

Catalog Number:
85730-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
243146G8
Reactivity:
human
GenBank:
BC024592
Immunogen Catalog Number:
Ag3859

Conjugate:
Unconjugated
Full name:
neurochondrin
Gene ID:
23154

Applications

Tested Applications:
Cytometric bead array

Range:
0.625-80 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP02072-2 targets Neurochondrin in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: NCDN Recombinant antibody, PBS Only (Capture) 85730-4-PBS (243146D9). 100 µg. Concentration 1 mg/mL.

Detection antibody: NCDN Recombinant antibody, PBS Only (Detector) 85730-2-PBS (243146G8). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP02072-2, Neurochondrin Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85730-4-PBS. Detection antibody: 85730-2-PBS. Standard: Ag3859. Range: 0.625-80 ng/mL