

# AKT3 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP02078-1

## Capture Antibody Information

**Catalog Number:**  
85769-1-PBS

**Host:**  
Rabbit

**Isotype:**  
IgG

**Purification Method:**  
Protein A purification

**Clone ID:**  
250046H5

**Reactivity:**  
human

**GenBank:**  
BC121154

**Immunogen Catalog Number:**  
Ag16298

**Conjugate:**  
Unconjugated

**Full name:**  
v-akt murine thymoma viral oncogene homolog 3 (protein kinase B, gamma)

**Gene ID:**  
10000

## Detection Antibody Information

**Catalog Number:**  
85769-2-PBS

**Host:**  
Rabbit

**Isotype:**  
IgG

**Purification Method:**  
Protein A purification

**Clone ID:**  
250046D12

**Reactivity:**  
human

**GenBank:**  
BC121154

**Immunogen Catalog Number:**  
Ag16298

**Conjugate:**  
Unconjugated

**Full name:**  
v-akt murine thymoma viral oncogene homolog 3 (protein kinase B, gamma)

**Gene ID:**  
10000

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.156-20 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP02078-1 targets AKT3 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: AKT3 Recombinant antibody, PBS Only (Capture) 85769-1-PBS (250046H5). 100 µg. Concentration 1 mg/mL.

Detection antibody: AKT3 Recombinant antibody, PBS Only (Detector) 85769-2-PBS (250046D12). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
**Storage buffer:**  
PBS only

# Selected Validation Data



Cytometric bead array standard curve of MP02078-1, AKT3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85769-1-PBS. Detection antibody: 85769-2-PBS. Standard: Ag16298. Range: 0.156-20 ng/mL