

LRG1 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP02084-1**

Capture Antibody Information

Catalog Number:

85762-2-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

243107G9

Reactivity:

human

GenBank:

NM_052972.3

Conjugate:

Unconjugated

Full name:

leucine-rich alpha-2-glycoprotein 1

Gene ID:

116844

Detection Antibody Information

Catalog Number:

85762-1-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

243107E11

Reactivity:

human

GenBank:

NM_052972.3

Conjugate:

Unconjugated

Full name:

leucine-rich alpha-2-glycoprotein 1

Gene ID:

116844

Applications

Tested Applications:

Cytometric bead array

Range:

3.125-200 ng/mL (Cytometric Bead Array)

Recommended Dilutions:

It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP02084-1 targets LRG1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: LRG1 Recombinant antibody, PBS Only (Capture) 85762-2-PBS (243107G9). 100 µg. Concentration 1 mg/mL.

Detection antibody: LRG1 Recombinant antibody, PBS Only (Detector) 85762-1-PBS (243107E11). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:


Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:

PBS only

Selected Validation Data



Cytometric bead array standard curve of MP02084-1, LRG1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85762-2-PBS. Detection antibody: 85762-1-PBS. Standard: Eg3083. Range: 3.125-200 ng/mL.