

OIP5 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP02091-2**

Capture Antibody Information

Catalog Number: 85757-4-PBS	Clone ID: 250073G5	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: Opa interacting protein 5
Isotype: IgG	GenBank: BC015050	Gene ID: 11339
Purification Method: Protein A purification	Immunogen Catalog Number: Ag2790	

Detection Antibody Information

Catalog Number: 85757-3-PBS	Clone ID: 250073E9	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: Opa interacting protein 5
Isotype: IgG	GenBank: BC015050	Gene ID: 11339
Purification Method: Protein A purification	Immunogen Catalog Number: Ag2790	

Applications

Tested Applications: Cytometric bead array	Range: 3.125-200 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
------------------------------------------------------	----------------------------------------------------------	-------------------------------------------------------------------------------------------------------------------------------------------

Product Information

MP02091-2 targets OIP5 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: OIP5 Recombinant antibody, PBS Only (Capture) 85757-4-PBS (250073G5). 100 µg. Concentration 1 mg/mL.

Detection antibody: OIP5 Recombinant antibody, PBS Only (Detector) 85757-3-PBS (250073E9). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP02091-2, OIP5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85757-4-PBS. Detection antibody: 85757-3-PBS. Standard: Ag2790. Range: 3.125-200 ng/mL.