

# TNFAIP8L2 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP02092-2**

## Capture Antibody Information

Catalog Number:  
85759-3-PBS  
Host:  
Rabbit  
Isotype:  
IgG  
Purification Method:  
Protein A purification

Clone ID:  
250047B6  
Reactivity:  
human  
GenBank:  
BC063014  
Immunogen Catalog Number:  
Ag7964

Conjugate:  
Unconjugated  
Full name:  
tumor necrosis factor, alpha-induced protein 8-like 2  
Gene ID:  
79626

## Detection Antibody Information

Catalog Number:  
85759-2-PBS  
Host:  
Rabbit  
Isotype:  
IgG  
Purification Method:  
Protein A purification

Clone ID:  
250047H4  
Reactivity:  
human  
GenBank:  
BC063014  
Immunogen Catalog Number:  
Ag7964

Conjugate:  
Unconjugated  
Full name:  
tumor necrosis factor, alpha-induced protein 8-like 2  
Gene ID:  
79626

## Applications

Tested Applications:  
Cytometric bead array

Range:  
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP02092-2 targets TNFAIP8L2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.  
Capture antibody: TNFAIP8L2 Recombinant antibody, PBS Only 85759-3-PBS (250047B6). 100 µg. Concentration 1 mg/mL.

Detection antibody: TNFAIP8L2 Recombinant antibody, PBS Only 85759-2-PBS (250047H4). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
Storage buffer:  
PBS only

# Selected Validation Data



Cytometric bead array standard curve of MP02092-2, TNFAIP8L2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85759-3-PBS. Detection antibody: 85759-2-PBS. Standard: Ag7964. Range: 0.781-100 ng/mL.