

# PDCD6 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP02102-1**

## Capture Antibody Information

Catalog Number: 85792-4-PBS	Clone ID: 250087C8	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: programmed cell death 6
Isotype: IgG	GenBank: BC012384	Gene ID: 10016
Purification Method: Protein A purification	Immunogen Catalog Number: Ag2949	

## Detection Antibody Information

Catalog Number: 85792-3-PBS	Clone ID: 250087C5	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: programmed cell death 6
Isotype: IgG	GenBank: BC012384	Gene ID: 10016
Purification Method: Protein A purification	Immunogen Catalog Number: Ag2949	

## Applications

Tested Applications: Cytometric bead array	Range: 0.781-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
---	---	--

## Product Information

MP02102-1 targets PDCD6 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: PDCD6 Recombinant antibody, PBS Only (Capture) 85792-4-PBS (250087C8). 100 µg. Concentration 1 mg/mL.

Detection antibody: PDCD6 Recombinant antibody, PBS Only (Detector) 85792-3-PBS (250087C5). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
Storage buffer:  
PBS only

# Selected Validation Data



Cytometric bead array standard curve of MP02102-1, PDCD6 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85792-4-PBS. Detection antibody: 85792-3-PBS. Standard: Ag2949. Range: 0.781-100 ng/mL.