

# Mouse IL-27 p28/IL-30 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP02119-1**

## Capture Antibody Information

**Catalog Number:**  
85795-3-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
242655D11  
**Reactivity:**  
mouse  
**GenBank:**  
**Immunogen Catalog Number:**  
Eg0897

**Conjugate:**  
Unconjugated  
**Full name:**  
interleukin 27  
**Gene ID:**  
246779

## Detection Antibody Information

**Catalog Number:**  
85795-2-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
242655D10  
**Reactivity:**  
mouse  
**GenBank:**

**Conjugate:**  
Unconjugated  
**Full name:**  
interleukin 27  
**Gene ID:**  
246779

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.781-50 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP02119-1 targets IL-27 p28/IL-30 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Mouse IL-27A p28 Recombinant antibody, PBS Only (Capture) 85795-3-PBS (242655D11). 100 µg. Concentration 1 mg/mL.

Detection antibody: Mouse IL-27A Recombinant antibody, PBS Only (Detector) 85795-2-PBS (242655D10). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
**Storage buffer:**  
PBS only

# Selected Validation Data



Cytometric bead array standard curve of MP02119-1, MOUSE IL-27 p28/IL-30 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85795-3-PBS. Detection antibody: 85795-2-PBS. Standard: Eg0897. Range: 0.781-50 ng/mL