

FAK Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP02124-2**

Capture Antibody Information

Catalog Number:
81334-4-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
242295E3
Reactivity:
human
GenBank:
BC028733
Immunogen Catalog Number:
Ag3331

Conjugate:
Unconjugated
Full name:
PTK2 protein tyrosine kinase 2
Gene ID:
5747

Detection Antibody Information

Catalog Number:
81334-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
242295B2
Reactivity:
human
GenBank:
BC028733
Immunogen Catalog Number:
Ag3331

Conjugate:
Unconjugated
Full name:
PTK2 protein tyrosine kinase 2
Gene ID:
5747

Applications

Tested Applications:
Cytometric bead array

Range:
0.156-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP02124-2 targets FAK in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: FAK Recombinant antibody, PBS Only (Capture) 81334-4-PBS (242295E3). 100 μ g. Concentration 1 mg/mL.

Detection antibody: FAK Recombinant antibody, PBS Only (Detector) 81334-3-PBS (242295B2). 100 μ g. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP02124-2, FAK Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 81334-4-PBS. Detection antibody: 81334-3-PBS. Standard: Ag3331. Range: 0.156-20 ng/mL