

Mouse uPAR/CD87 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP02130-1**

Capture Antibody Information

| | | |
|---|--|--|
| Catalog Number: 85824-4-PBS | Clone ID: 242976C11 | Conjugate: Unconjugated |
| Host: Rabbit | Reactivity: mouse | Full name: plasminogen activator, urokinase receptor |
| Isotype: IgG | GenBank: NM_011113.4 | Gene ID: 18793 |
| Purification Method: Protein A purification | Immunogen Catalog Number: Eg3156 | |

Detection Antibody Information

| | | |
|---|--|--|
| Catalog Number: 85824-3-PBS | Clone ID: 242976C2 | Conjugate: Unconjugated |
| Host: Rabbit | Reactivity: mouse | Full name: plasminogen activator, urokinase receptor |
| Isotype: IgG | GenBank: NM_011113.4 | Gene ID: 18793 |
| Purification Method: Protein A purification | Immunogen Catalog Number: Eg3156 | |

Applications

| | | |
|--|---|---|
| Tested Applications: Cytometric bead array | Range: 0.625-80 ng/mL (Cytometric Bead Array) | Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results. |
|--|---|---|

Product Information

MP02130-1 targets uPAR/CD87 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Mouse uPAR/CD87 Recombinant antibody, PBS Only 85824-4-PBS (242976C11). 100 µg. Concentration 1 mg/mL.

Detection antibody: Mouse uPAR/CD87 Recombinant antibody, PBS Only 85824-3-PBS (242976C2). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP02130-1, MOUSE uPAR/CD87 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85824-4-PBS. Detection antibody: 85824-3-PBS. Standard: Eg3156. Range: 0.625-80 ng/mL