

For Research Use Only

Mouse CCL19/MIP-3 beta Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP02165-2

Capture Antibody Information

Catalog Number:
85861-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
242867E7
Reactivity:
mouse
GenBank:
NM_011888.2
Immunogen Catalog Number:
Eg3124

Conjugate:
Unconjugated
Full name:
chemokine (C-C motif) ligand 19
Gene ID:
24047

Detection Antibody Information

Catalog Number:
85861-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
242867C12
Reactivity:
mouse
GenBank:
NM_011888.2
Immunogen Catalog Number:
Eg3124

Conjugate:
Unconjugated
Full name:
chemokine (C-C motif) ligand 19
Gene ID:
24047

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP02165-2 targets CCL19/MIP-3 beta in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Mouse Ccl19 Recombinant antibody, PBS Only (Capture) 85861-3-PBS (242867E7). 100 µg. Concentration 1 mg/mL.

Detection antibody: Mouse Ccl19 Recombinant antibody, PBS Only 85861-2-PBS (242867C12). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptgcn.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP02165-2, MOUSE CCL19/MIP-3 beta Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85861-3-PBS. Detection antibody: 85861-2-PBS. Standard: Eg3124. Range: 0.781-100 ng/mL.