

Mouse CD62L Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP02196-1

Capture Antibody Information

Catalog Number: 85903-3-PBS	Clone ID: 250271A4	Conjugate: Unconjugated
Host: Rabbit	Reactivity: mouse	Full name: selectin, lymphocyte
Isotype: IgG	GenBank: NM_011346.2	Gene ID: 20343
Purification Method: Protein A purification	Immunogen Catalog Number: Eg3845	

Detection Antibody Information

Catalog Number: 85903-2-PBS	Clone ID: 250271H3	Conjugate: Unconjugated
Host: Rabbit	Reactivity: mouse	Full name: selectin, lymphocyte
Isotype: IgG	GenBank: NM_011346.2	Gene ID: 20343
Purification Method: Protein A purification	Immunogen Catalog Number: Eg3845	

Applications

Tested Applications: Cytometric bead array	Range: 0.781-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
------------------------------------------------------	----------------------------------------------------------	-------------------------------------------------------------------------------------------------------------------------------------------

Product Information

MP02196-1 targets CD62L in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Mouse CD62L Recombinant antibody, PBS Only (Capture) 85903-3-PBS (250271A4). 100 µg. Concentration 1 mg/mL.

Detection antibody: Mouse CD62L Recombinant antibody, PBS Only (Detector) 85903-2-PBS (250271H3). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP02196-1, MOUSE CD62L Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85903-3-PBS. Detection antibody: 85903-2-PBS. Standard: Eg3845. Range: 0.781-100 ng/mL.