

Mouse Ccl8 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP02245-1**

Capture Antibody Information

Catalog Number: 86081-2-PBS	Clone ID: 250174H8	Conjugate: Unconjugated
Host: Rabbit	Reactivity: mouse	Full name: chemokine (C-C motif) ligand 8
Isotype: IgG	GenBank: NM_021443.3	Gene ID: 20307
Purification Method: Protein A purification	Immunogen Catalog Number: Eg2771	

Detection Antibody Information

Catalog Number: 86081-1-PBS	Clone ID: 250174H6	Conjugate: Unconjugated
Host: Rabbit	Reactivity: mouse	Full name: chemokine (C-C motif) ligand 8
Isotype: IgG	GenBank: NM_021443.3	Gene ID: 20307
Purification Method: Protein A purification	Immunogen Catalog Number: Eg2771	

Applications

Tested Applications: Cytometric bead array	Range: 0.156-20 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	---	---

Product Information

MP02245-1 targets Ccl8 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Mouse Ccl8 Recombinant antibody, PBS Only (Capture) 86081-2-PBS (250174H8). 100 µg. Concentration 1 mg/mL.

Detection antibody: Mouse Ccl8 Recombinant antibody, PBS Only (Detector) 86081-1-PBS (250174H6). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



The mean Ccl8 concentration was determined to be 1317.1 ng/mL in mouse serum.



Cytometric bead array standard curve of MP02245-1, MOUSE Ccl8 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86081-2-PBS. Detection antibody: 86081-1-PBS. Standard: Eg2771. Range: 0.156-20 ng/mL.