

Neuropilin 2 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP02281-1**

Capture Antibody Information

Catalog Number:
86126-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
250688B8
Reactivity:
human
GenBank:
NM_201266.1
Immunogen Catalog Number:
Eg3125

Conjugate:
Unconjugated
Full name:
neuropilin 2
Gene ID:
8828

Detection Antibody Information

Catalog Number:
86126-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
250688E10
Reactivity:
human
GenBank:
NM_201266.1
Immunogen Catalog Number:
Eg3125

Conjugate:
Unconjugated
Full name:
neuropilin 2
Gene ID:
8828

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP02281-1 targets Neuropilin 2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: NRP2 Recombinant antibody, PBS Only (Capture) 86126-2-PBS (250688B8). 100 µg. Concentration 1 mg/mL.

Detection antibody: NRP2 Recombinant antibody, PBS Only (Detector) 86126-1-PBS (250688E10). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



The mean Neuropilin 2 concentration was determined to be 15.7 ng/mL in human serum.



Cytometric bead array standard curve of MP02281-1, Neuropilin 2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86126-2-PBS. Detection antibody: 86126-1-PBS. Standard: Eg3125. Range: 0.781-100 ng/mL