

# LPP Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP02282-2**

## Capture Antibody Information

**Catalog Number:**  
86117-2-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
250773F6  
**Reactivity:**  
human  
**GenBank:**  
BC130584  
**Immunogen Catalog Number:**  
Ag19205

**Conjugate:**  
Unconjugated  
**Full name:**  
LIM domain containing preferred translocation partner in lipoma  
**Gene ID:**  
4026

## Detection Antibody Information

**Catalog Number:**  
86117-4-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
250773B3  
**Reactivity:**  
human  
**GenBank:**  
BC130584  
**Immunogen Catalog Number:**  
Ag19205

**Conjugate:**  
Unconjugated  
**Full name:**  
LIM domain containing preferred translocation partner in lipoma  
**Gene ID:**  
4026

## Applications

**Tested Applications:**  
Sandwich ELISA

**Range:**  
78.1-5000 pg/mL (Sandwich ELISA)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP02282-2 targets LPP in immunoassays as a matched antibody pair. Validated in Sandwich ELISA.

Capture antibody: LPP Recombinant antibody, PBS Only (Capture) 86117-2-PBS (250773F6). 100 µg. Concentration 1 mg/mL.

Detection antibody: LPP Recombinant antibody, PBS Only (Detector) 86117-4-PBS (250773B3). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
**Storage buffer:**  
PBS only

## Selected Validation Data



Sandwich ELISA standard curve of MP02282-2, Human LPP Recombinant Matched Antibody Pair - PBS only. 86117-2-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag19205. 86117-4-PBS was HRP conjugated as the detection antibody. Range: 78.1-5000 pg/mL.



The mean LPP concentration was determined to be 43,872.1 pg/mL in MCF-7 cell extract based on a 1.2 mg/mL extract load and 4,111.8 pg/mL in NIH/3T3 cell extract based on a 1.4 mg/mL extract load.