

# Progranulin/PGRN Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP02283-1**

## Capture Antibody Information

Catalog Number: 86118-2-PBS	Clone ID: 250681G2	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: granulin
Isotype: IgG	GenBank: NM_002087.3	Gene ID: 2896
Purification Method: Protein A purification	Immunogen Catalog Number: Eg2954	

## Detection Antibody Information

Catalog Number: 86118-1-PBS	Clone ID: 250681G12	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: granulin
Isotype: IgG	GenBank: NM_002087.3	Gene ID: 2896
Purification Method: Protein A purification	Immunogen Catalog Number: Eg2954	

## Applications

Tested Applications: Cytometric bead array	Range: 0.625-80 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
---	--	--

## Product Information

MP02283-1 targets Progranulin/PGRN in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: GRN Recombinant antibody, PBS Only (Capture) 86118-2-PBS (250681G2). 100 µg. Concentration 1 mg/mL.

Detection antibody: GRN Recombinant antibody, PBS Only (Detector) 86118-1-PBS (250681G12). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

## Storage


Storage:  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage buffer:  
PBS only

## Selected Validation Data



The mean Progranulin/PGRN concentration was determined to be 64.7 ng/mL in human serum.



Cytometric bead array standard curve of MP02283-1, Progranulin/PGRN Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86118-2-PBS. Detection antibody: 86118-1-PBS. Standard: Eg2954. Range: 0.625-80 ng/mL.