

# Frataxin Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP02292-1**

## Capture Antibody Information

Catalog Number:  
86195-2-PBS  
Host:  
Rabbit  
Isotype:  
IgG  
Purification Method:  
Protein A purification

Clone ID:  
250761C1  
Reactivity:  
human  
GenBank:  
NM\_000144.5  
Immunogen Catalog Number:  
Eg4622

Conjugate:  
Unconjugated  
Full name:  
frataxin  
Gene ID:  
2395

## Detection Antibody Information

Catalog Number:  
86195-1-PBS  
Host:  
Rabbit  
Isotype:  
IgG  
Purification Method:  
Protein A purification

Clone ID:  
250761A7  
Reactivity:  
human  
GenBank:  
NM\_000144.5  
Immunogen Catalog Number:  
Eg4622

Conjugate:  
Unconjugated  
Full name:  
frataxin  
Gene ID:  
2395

## Applications

Tested Applications:  
Cytometric bead array

Range:  
0.313-40 ng/mL (Cytometric Bead Array)

Recommended Dilutions:  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP02292-1 targets Frataxin in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Frataxin Recombinant antibody, PBS Only (Capture) 86195-2-PBS (250761C1). 100 µg. Concentration 1 mg/mL.

Detection antibody: Frataxin Recombinant antibody, PBS Only (Detector) 86195-1-PBS (250761A7). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


## Storage

Storage:  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
Storage buffer:  
PBS only

## Selected Validation Data



The mean Frataxin concentration was determined to be 779.5 ng/mL in HepG2 cell extract based on a 0.9 mg/mL extract load.



Cytometric bead array standard curve of MP02292-1, Frataxin Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86195-2-PBS. Detection antibody: 86195-1-PBS. Standard: Eg4622. Range: 0.313-40 ng/mL.