

Cytokeratin 19 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP02308-1**

Capture Antibody Information

Catalog Number:

86250-3-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

250894D12

Reactivity:

human

GenBank:

BC002539

Conjugate:

Unconjugated

Full name:

keratin 19

Gene ID:

3880

Detection Antibody Information

Catalog Number:

86250-2-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

250894E8

Reactivity:

human

GenBank:

BC002539

Conjugate:

Unconjugated

Full name:

keratin 19

Gene ID:

3880

Applications

Tested Applications:

Cytometric bead array

Range:

1.56-200 ng/mL (Cytometric Bead Array)

Recommended Dilutions:

It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP02308-1 targets Cytokeratin 19 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: KRT19 Recombinant antibody, PBS Only (Capture) 86250-3-PBS (250894D12). 100 µg. Concentration 1 mg/mL.

Detection antibody: KRT19 Recombinant antibody, PBS Only (Detector) 86250-2-PBS (250894E8). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:


Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C


Storage buffer:

PBS only

Selected Validation Data



The mean Cytokeratin 19 concentration was determined to be 183.9 ng/mL in MCF-7 cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP02308-1, Cytokeratin 19 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86250-3-PBS. Detection antibody: 86250-2-PBS. Standard: Ag7407. Range: 1.56-200 ng/mL.