

# FKBP1A Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP02313-1**

## Capture Antibody Information

<b>Catalog Number:</b> 86249-3-PBS	<b>Clone ID:</b> 250923D6	<b>Conjugate:</b> Unconjugated
<b>Host:</b> Rabbit	<b>Reactivity:</b> human	<b>Full name:</b> FK506 binding protein 1A, 12kDa
<b>Isotype:</b> IgG	<b>GenBank:</b> BC001925	<b>Gene ID:</b> 2280
<b>Purification Method:</b> Protein A purification	<b>Immunogen Catalog Number:</b> Ag0406	

## Detection Antibody Information

<b>Catalog Number:</b> 86249-2-PBS	<b>Clone ID:</b> 250923G1	<b>Conjugate:</b> Unconjugated
<b>Host:</b> Rabbit	<b>Reactivity:</b> human	<b>Full name:</b> FK506 binding protein 1A, 12kDa
<b>Isotype:</b> IgG	<b>GenBank:</b> BC001925	<b>Gene ID:</b> 2280
<b>Purification Method:</b> Protein A purification	<b>Immunogen Catalog Number:</b> Ag0406	

## Applications

<b>Tested Applications:</b> Cytometric bead array	<b>Range:</b> 0.781-100 ng/mL (Cytometric Bead Array)	<b>Recommended Dilutions:</b> It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
------------------------------------------------------	----------------------------------------------------------	-------------------------------------------------------------------------------------------------------------------------------------------

## Product Information

MP02313-1 targets FKBP1A in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: FKBP1A Recombinant antibody, PBS Only (Capture) 86249-3-PBS (250923D6). 100 µg. Concentration 1 mg/mL.

Detection antibody: FKBP1A Recombinant antibody, PBS Only (Detector) 86249-2-PBS (250923G1). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
**Storage buffer:**  
PBS only

## Selected Validation Data



The mean FKBP1A concentration was determined to be 149.2 ng/mL in Jurkat cell extract based on a 2.2 mg/mL extract load.



Cytometric bead array standard curve of MP02313-1, FKBP1A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86249-3-PBS. Detection antibody: 86249-2-PBS. Standard: Ag0406. Range: 0.781-100 ng/mL.