

For Research Use Only

BCORL1 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP02366-1

Capture Antibody Information

Catalog Number:
86290-1-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
251013E2

Reactivity:
human

GenBank:
BC156742

Immunogen Catalog Number:
Ag19522

Conjugate:
Unconjugated

Full name:
BCL6 co-repressor-like 1

Gene ID:
63035

Detection Antibody Information

Catalog Number:
86290-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
251013G12

Reactivity:
human

GenBank:
BC156742

Immunogen Catalog Number:
Ag19522

Conjugate:
Unconjugated

Full name:
BCL6 co-repressor-like 1

Gene ID:
63035

Applications

Tested Applications:
Cytometric bead array

Range:
0.117-30 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP02366-1 targets BCORL1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: BCORL1 Recombinant monoclonal antibody, PBS Only (Capture) 86290-1-PBS (251013E2). 100 μ g. Concentration 1 mg/mL.

Detection antibody: BCORL1 Recombinant monoclonal antibody, PBS Only (Detector) 86290-2-PBS (251013G12). 100 μ g. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

For technical support and original validation data for this product please contact:



T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MPO2366-1, BCORL1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86290-1-PBS. Detection antibody: 86290-2-PBS. Standard: Ag19522. Range: 0.117-30 ng/mL

The mean BCORL1 concentration was determined to be 23.1 ng/mL in K-562 cell extract based on a 3.0 mg/mL extract load.