

For Research Use Only



HOXA10 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP02388-1

Capture Antibody Information

Catalog Number:
85951-4-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
250418C10
Reactivity:
human
GenBank:
BC013971
Immunogen Catalog Number:
Ag24156

Conjugate:
Unconjugated
Full name:
homeobox A10
Gene ID:
3206

Detection Antibody Information

Catalog Number:
85951-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
250418C7
Reactivity:
human
GenBank:
BC013971
Immunogen Catalog Number:
Ag24156

Conjugate:
Unconjugated
Full name:
homeobox A10
Gene ID:
3206

Applications

Tested Applications:
Sandwich ELISA

Range:
39.1-2500 pg/mL (Sandwich ELISA)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP02388-1 targets HOXA10 in immunoassays as a matched antibody pair. Validated in Sandwich ELISA.

Capture antibody: HOXA10 Recombinant monoclonal antibody, PBS Only (Capture) 85951-4-PBS (250418C10). 100 µg. Concentration 1 mg/mL.

Detection antibody: HOXA10 Recombinant monoclonal antibody, PBS Only (Detector) 85951-3-PBS (250418C7). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP02388-1, Human HOXA10 Recombinant Matched Antibody Pair - PBS only. 85951-4-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag24156. 85951-3-PBS was HRP conjugated as the detection antibody. Range: 39.1-2500 pg/mL.



The mean HOXA10 concentration was determined to be 21.8 ng/mL in K-562 cell extract based on a 3.0 mg/mL extract load and 56.2 ng/mL in HEK-293T cell extract based on a 1.3 mg/mL extract load.