

REEP4 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP02417-1

Capture Antibody Information

Catalog Number: 86164-5-PBS	Clone ID: 250770D9	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: receptor accessory protein 4
Isotype: IgG	GenBank: BC013048	Gene ID: 80346
Purification Method: Protein A purification	Immunogen Catalog Number: Ag24569	

Detection Antibody Information

Catalog Number: 86164-4-PBS	Clone ID: 250770A8	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: receptor accessory protein 4
Isotype: IgG	GenBank: BC013048	Gene ID: 80346
Purification Method: Protein A purification	Immunogen Catalog Number: Ag24569	

Applications

Tested Applications: Cytometric bead array	Range: 0.781-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

Product Information

MP02417-1 targets REEP4 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: REEP4 Recombinant antibody, PBS Only (Capture) 86164-5-PBS (250770D9). 100 µg. Concentration 1 mg/mL.

Detection antibody: REEP4 Recombinant antibody, PBS Only (Detector) 86164-4-PBS (250770A8). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



The mean REEP4 concentration was determined to be 5.4 ng/mL in HeLa cell extract based on a 1.5 mg/mL extract load.



Cytometric bead array standard curve of MP02417-1, REEP4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86164-5-PBS. Detection antibody: 86164-4-PBS. Standard: Ag24569. Range: 0.781-100 ng/mL.