

For Research Use Only

ARD1A Recombinant Matched Antibody Pair, PBS Only



Catalog Number: MP02442-1

Capture Antibody Information

Catalog Number:
86450-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
251237B11
Reactivity:
human
GenBank:
BC000308
Immunogen Catalog Number:
Ag6519

Conjugate:
Unconjugated
Full name:
ARD1 homolog A, N-acetyltransferase (S. cerevisiae)
Gene ID:
8260

Detection Antibody Information

Catalog Number:
86450-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
251237D2
Reactivity:
human
GenBank:
BC000308
Immunogen Catalog Number:
Ag6519

Conjugate:
Unconjugated
Full name:
ARD1 homolog A, N-acetyltransferase (S. cerevisiae)
Gene ID:
8260

Applications

Tested Applications:
Cytometric bead array

Range:
0.469-60 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP02442-1 targets ARD1A in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: ARD1A Recombinant antibody, PBS Only (Capture) 86450-2-PBS (251237B11). 100 µg. Concentration 1 mg/mL.

Detection antibody: ARD1A Recombinant antibody, PBS Only (Detector) 86450-1-PBS (251237D2). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean ARD1A concentration was determined to be 13.8 ng/mL in MCF-7 cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP02442-1, ARD1A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86450-2-PBS. Detection antibody: 86450-1-PBS. Standard: Ag6519. Range: 0.469-60 ng/mL.