

GSPT2 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP02477-2**

Capture Antibody Information

Catalog Number: 86508-1-PBS	Clone ID: 251224C2	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: G1 to S phase transition 2
Isotype: IgG	GenBank: BC036077	Gene ID: 23708
Purification Method: Protein A purification	Immunogen Catalog Number: Ag3641	

Detection Antibody Information

Catalog Number: 86508-4-PBS	Clone ID: 251224G9	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: G1 to S phase transition 2
Isotype: IgG	GenBank: BC036077	Gene ID: 23708
Purification Method: Protein A purification	Immunogen Catalog Number: Ag3641	

Applications

Tested Applications: Sandwich ELISA	Range: 0.156-10 ng/mL (Sandwich ELISA)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	---	--

Product Information

MP02477-2 targets GSPT2 in immunoassays as a matched antibody pair. Validated in Sandwich ELISA.

Capture antibody: GSPT2 Recombinant monoclonal antibody, PBS Only (Capture/Detector) 86508-1-PBS (251224C2). 100 µg. Concentration 1 mg/mL.

Detection antibody: GSPT2 Recombinant monoclonal antibody, PBS Only (Detector) 86508-4-PBS (251224G9). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Sandwich ELISA standard curve of MP02477-2, Human GSPT2 Recombinant Matched Antibody Pair - PBS only. 86508-1-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag3641. 86508-4-PBS was HRP conjugated as the detection antibody. Range: 0.156-10 ng/mL.



The mean GSPT2 concentration was determined to be 122.43 ng/mL in MCF-7 cell extract based on a 1.20 mg/mL extract load.