

MSRA Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP02504-1**

Capture Antibody Information

Catalog Number: 86358-4-PBS	Clone ID: 251025C2	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: methionine sulfoxide reductase A
Isotype: IgG	GenBank: BC054033	Gene ID: 4482
Purification Method: Protein A purification	Immunogen Catalog Number: Ag6053	

Detection Antibody Information

Catalog Number: 86358-3-PBS	Clone ID: 251025F1	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: methionine sulfoxide reductase A
Isotype: IgG	GenBank: BC054033	Gene ID: 4482
Purification Method: Protein A purification	Immunogen Catalog Number: Ag6053	

Applications

Tested Applications: Sandwich ELISA	Range: 0.156–10 ng/mL (Sandwich ELISA)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
---	--	---

Product Information

MP02504-1 targets MSRA in immunoassays as a matched antibody pair. Validated in Sandwich ELISA.

Capture antibody: MSRA Recombinant monoclonal antibody, PBS Only (Capture) 86358-4-PBS (251025C2). 100 µg. Concentration 1 mg/mL.

Detection antibody: MSRA Recombinant monoclonal antibody, PBS Only (Detector) 86358-3-PBS (251025F1). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Sandwich ELISA standard curve of MP02504-1, Human MSRA Recombinant Matched Antibody Pair - PBS only. 86358-4-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag6053. 86358-3-PBS was HRP conjugated as the detection antibody. Range: 0.156-10 ng/mL



The mean MSRA concentration was determined to be 1.82 ng/mL in HT-1080 cell extract based on a 1.4 mg/mL extract load.