

Syntaxin 8 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP02538-1**

Capture Antibody Information

Catalog Number:
86580-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
251516H4
Reactivity:
human
GenBank:
BC009713
Immunogen Catalog Number:
Ag2850

Conjugate:
Unconjugated
Full name:
syntaxin 8
Gene ID:
9482

Detection Antibody Information

Catalog Number:
86580-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
251516G10
Reactivity:
human
GenBank:
BC009713
Immunogen Catalog Number:
Ag2850

Conjugate:
Unconjugated
Full name:
syntaxin 8
Gene ID:
9482

Applications

Tested Applications:
Cytometric bead array

Range:
0.469-60 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP02538-1 targets Syntaxin 8 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: STX8 Recombinant antibody, PBS Only (Capture) 86580-2-PBS (251516H4). 100 µg. Concentration 1 mg/mL.

Detection antibody: STX8 Recombinant antibody, PBS Only (Detector) 86580-1-PBS (251516G10). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



The mean Syntaxin 8 concentration was determined to be 28.1 ng/mL in HepG2 cell extract based on a 1.2 mg/mL extract load, 43.4 ng/mL in A375 cell extract based on a 1.6 mg/mL extract load.



Cytometric bead array standard curve of MP02538-1, Syntaxin 8 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86580-2-PBS. Detection antibody: 86580-1-PBS. Standard: Ag2850. Range: 0.469-60 ng/mL.