

Ataxin 2 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP02584-1

Capture Antibody Information

Catalog Number: 80919-2-PBS	Clone ID: 250888E4	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human, mouse	Full name: ataxin 2
Isotype: IgG	GenBank: BC114546	Gene ID: 6311
Purification Method: Protein A purification	Immunogen Catalog Number: Ag16470	

Detection Antibody Information

Catalog Number: 80919-4-PBS	Clone ID: 250888A7	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: ataxin 2
Isotype: IgG	GenBank: BC114546	Gene ID: 6311
Purification Method: Protein A purification	Immunogen Catalog Number: Ag16470	

Applications

Tested Applications: Sandwich ELISA	Range: 62.5-4000 pg/mL (Sandwich ELISA)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
---	---	---

Product Information

MP02584-1 targets Ataxin 2 in immunoassays as a matched antibody pair. Validated in Sandwich ELISA.

Capture antibody: Ataxin 2 Recombinant antibody, PBS Only (Capture) 80919-2-PBS (250888E4). 100 µg. Concentration 1 mg/mL.

Detection antibody: Ataxin 2 Recombinant antibody, PBS Only (Detector) 80919-4-PBS (250888A7). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Sandwich ELISA standard curve of MP02584-1, Human Ataxin 2 Recombinant Matched Antibody Pair - PBS only. 80919-2-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag16470. 80919-4-PBS was HRP conjugated as the detection antibody. Range: 62.5-4000 pg/mL



The mean Ataxin 2 concentration was determined to be 2.25 ng/mL in HepG2 cell extract based on a 1.3 mg/mL extract load.