

For Research Use Only

MBD3 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP02598-1

Capture Antibody Information

Catalog Number:
86681-3-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
250922A5

Reactivity:
human

GenBank:
BC043619

Immunogen Catalog Number:
Ag5531

Conjugate:
Unconjugated

Full name:
methyl-CpG binding domain protein 3

Gene ID:
53615

Detection Antibody Information

Catalog Number:
86681-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
250922E11

Reactivity:
human

GenBank:
BC043619

Immunogen Catalog Number:
Ag5531

Conjugate:
Unconjugated

Full name:
methyl-CpG binding domain protein 3

Gene ID:
53615

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP02598-1 targets MBD3 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: MBD3 Recombinant monoclonal antibody, PBS Only (Capture) 86681-3-PBS (250922A5). 100 μ g. Concentration 1 mg/mL.

Detection antibody: MBD3 Recombinant monoclonal antibody, PBS Only (Detector) 86681-2-PBS (250922E11). 100 μ g. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:

PBS only

For technical support and original validation data for this product please contact:



T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP02598-1, MBD3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86681-3-PBS. Detection antibody: 86681-2-PBS. Standard: Ag5531. Range: 0.781-100 ng/mL.

The mean MBD3 concentration was determined to be 19.2 ng/mL in Y79 cell extract based on a 1.1 mg/mL extract load, 28.9 ng/mL in HeLa cell extract based on a 1.2 mg/mL extract load, 89.6 ng/mL in HEK-293 cell extract based on a 3.2 mg/mL extract load.