

For Research Use Only

Mouse PTPRZ1 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP02599-1

Capture Antibody Information

Catalog Number:
86687-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
251663D6
Reactivity:
mouse
GenBank:
NM_001081306.1
Immunogen Catalog Number:
Eg6825

Conjugate:
Unconjugated
Full name:
protein tyrosine phosphatase, receptor type Z, polypeptide 1
Gene ID:
19283

Detection Antibody Information

Catalog Number:
86687-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
251663C5
Reactivity:
mouse
GenBank:
NM_001081306.1
Immunogen Catalog Number:
Eg6825

Conjugate:
Unconjugated
Full name:
protein tyrosine phosphatase, receptor type Z, polypeptide 1
Gene ID:
19283

Applications

Tested Applications:
Cytometric bead array

Range:
0.312-40 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP02599-1 targets PTPRZ1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Mouse PTPRZ Recombinant monoclonal antibody, PBS Only (Capture/Detector) 86687-3-PBS (251663D6). 100 μ g. Concentration 1 mg/mL.

Detection antibody: Mouse PTPRZ Recombinant monoclonal antibody, PBS Only (Detector) 86687-2-PBS (251663C5). 100 μ g. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:



T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP02599-1, MOUSE PTPRZ1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86687-3-PBS. Detection antibody: 86687-2-PBS. Standard: Eg6825. Range: 0.312-40 ng/mL.

The mean PTPRZ1 concentration was determined to be 7.3 ng/mL in mouse cerebellum tissue extract based on a 2.9 mg/mL extract load.