

Mouse Cd2ap Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP02649-1**

Capture Antibody Information

Catalog Number: 86727-3-PBS	Clone ID: 251428H6	Conjugate: Unconjugated
Host: Rabbit	Reactivity: mouse	Full name: CD2-associated protein
Isotype: IgG	GenBank: BC157901	Gene ID: 12488
Purification Method: Protein A purification	Immunogen Catalog Number: Ag0622	

Detection Antibody Information

Catalog Number: 86727-2-PBS	Clone ID: 251428G6	Conjugate: Unconjugated
Host: Rabbit	Reactivity: mouse	Full name: CD2-associated protein
Isotype: IgG	GenBank: BC157901	Gene ID: 12488
Purification Method: Protein A purification	Immunogen Catalog Number: Ag0622	

Applications

Tested Applications: Cytometric bead array	Range: 0.781-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

Product Information

MP02649-1 targets Cd2ap in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Mouse Cd2ap Recombinant antibody, PBS Only (Capture) 86727-3-PBS (251428H6). 100 µg. Concentration 1 mg/mL.

Detection antibody: Mouse Cd2ap Recombinant antibody, PBS Only (Detector) 86727-2-PBS (251428G6). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



The mean Cd2ap concentration was determined to be 948.7 ng/mL in Mouse kidney tissue extract based on a 8.4 mg/mL extract load.



Cytometric bead array standard curve of MP02649-1, MOUSE Cd2ap Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86727-3-PBS. Detection antibody: 86727-2-PBS. Standard: Ag0622. Range: 0.781-100 ng/mL.