

For Research Use Only

Rat NCAM1/CD56 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP02663-2

Capture Antibody Information

Catalog Number:
86770-4-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
250717F7
Reactivity:
rat
GenBank:
NM_031521.1
Immunogen Catalog Number:
Eg3243

Conjugate:
Unconjugated
Full name:
neural cell adhesion molecule 1
Gene ID:
24586

Detection Antibody Information

Catalog Number:
86770-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
250717F6
Reactivity:
rat
GenBank:
NM_031521.1
Immunogen Catalog Number:
Eg3243

Conjugate:
Unconjugated
Full name:
neural cell adhesion molecule 1
Gene ID:
24586

Applications

Tested Applications:
Sandwich ELISA

Range:
0.781-50 ng/mL (Sandwich ELISA)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP02663-2 targets NCAM1/CD56 in immunoassays as a matched antibody pair. Validated in Sandwich ELISA.

Capture antibody: Rat NCAM-1/CD56 Recombinant monoclonal antibody, PBS Only (Capture) 86770-4-PBS (250717F7). 100 µg. Concentration 1 mg/mL.

Detection antibody: Rat NCAM-1/CD56 Recombinant monoclonal antibody, PBS Only (Detector) 86770-3-PBS (250717F6). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP02663-2, Rat NCAM1/CD56 Recombinant Matched Antibody Pair - PBS only. 86770-4-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg3243. 86770-3-PBS was HRP conjugated as the detection antibody. Range: 0.781-50 ng/mL.



The mean NCAM1/CD56 concentration was determined to be 89.12 ng/mL in C6 cell extract based on a 1.20 mg/mL extract load and 3,090.35 ng/mL in rat brain tissue extract based on a 8.2 mg/mL extract load.