

For Research Use Only

RGS7 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP02701-1

Capture Antibody Information

Catalog Number:
86846-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
251849E9

Reactivity:
human

GenBank:
BC022009

Immunogen Catalog Number:
Ag4982

Conjugate:
Unconjugated

Full name:
regulator of G-protein signaling 7

Gene ID:
6000

Detection Antibody Information

Catalog Number:
86846-1-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
251849C4

Reactivity:
human

GenBank:
BC022009

Immunogen Catalog Number:
Ag4982

Conjugate:
Unconjugated

Full name:
regulator of G-protein signaling 7

Gene ID:
6000

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP02701-1 targets RGS7 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: RGS7 Recombinant monoclonal antibody, PBS Only (Capture) 86846-2-PBS (251849E9). 100 μ g. Concentration 1 mg/mL.

Detection antibody: RGS7 Recombinant monoclonal antibody, PBS Only (Detector) 86846-1-PBS (251849C4). 100 μ g. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP02701-1, RGS7 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86846-2-PBS. Detection antibody: 86846-1-PBS. Standard: Ag4982. Range: 0.781-100 ng/mL.

The mean RGS7 concentration was determined to be 5.9 ng/mL in SH-SY5Y cell extract based on a 1.2 mg/mL extract load, 237.6 ng/mL in Mouse brain tissue extract based on a 3.2 mg/mL extract load.