

HDAC2 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP02705-1

Capture Antibody Information

Catalog Number: 81497-3-PBS	Clone ID: 251034C1	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: histone deacetylase 2
Isotype: IgG	GenBank: BC031055	Gene ID: 3066
Purification Method: Protein A purification	Immunogen Catalog Number: Ag3607	

Detection Antibody Information

Catalog Number: 81497-2-PBS	Clone ID: 251034E8	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: histone deacetylase 2
Isotype: IgG	GenBank: BC031055	Gene ID: 3066
Purification Method: Protein A purification	Immunogen Catalog Number: Ag3607	

Applications

Tested Applications: Sandwich ELISA	Range: 0.625-40 ng/mL (Sandwich ELISA)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	---	--

Product Information

MP02705-1 targets HDAC2 in immunoassays as a matched antibody pair. Validated in Sandwich ELISA.

Capture antibody: HDAC2 Recombinant monoclonal antibody, PBS Only (Capture) 81497-3-PBS (251034C1). 100 µg. Concentration 1 mg/mL.

Detection antibody: HDAC2 Recombinant monoclonal antibody, PBS Only (Detector) 81497-2-PBS (251034E8). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Sandwich ELISA standard curve of MP02705-1, Human HDAC2 Recombinant Matched Antibody Pair - PBS only. 81497-3-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag3607. 81497-2-PBS was HRP conjugated as the detection antibody. Range: 0.625-40 ng/mL



The mean HDAC2 concentration was determined to be 62.00 ng/mL in HEK-293 cell extract based on a 1.2 mg/mL extract load