

For Research Use Only

CPEB3 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP02708-1

Capture Antibody Information

Catalog Number:
86870-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
251275E1

Reactivity:
human

GenBank:
BC036444

Immunogen Catalog Number:
Ag3351

Conjugate:
Unconjugated

Full name:
cytoplasmic polyadenylation
element binding protein 3

Gene ID:
22849

Detection Antibody Information

Catalog Number:
86870-1-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
251275B3

Reactivity:
human

GenBank:
BC036444

Immunogen Catalog Number:
Ag3351

Conjugate:
Unconjugated

Full name:
cytoplasmic polyadenylation
element binding protein 3

Gene ID:
22849

Applications

Tested Applications:
Sandwich ELISA

Range:
39.1-2500 pg/mL (Sandwich ELISA)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP02708-1 targets CPEB3 in immunoassays as a matched antibody pair. Validated in Sandwich ELISA.

Capture antibody: CPEB3 Recombinant monoclonal antibody, PBS Only (Capture) 86870-2-PBS (251275E1). 100 µg. Concentration 1 mg/mL.

Detection antibody: CPEB3 Recombinant monoclonal antibody, PBS Only (Detector) 86870-1-PBS (251275B3). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP02708-1, Human CPEB3 Recombinant Matched Antibody Pair - PBS only. 86870-2-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag3351. 86870-1-PBS was HRP conjugated as the detection antibody. Range: 39.1-2500 pg/mL



The mean CPEB3 concentration was determined to be 7.00 ng/mL in HepG2 cell extract based on a 1.3 mg/mL extract load and 8.20 ng/mL in HuH-7 cell extract based on a 2.1 mg/mL extract load.