

WIBG Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP02733-1

Capture Antibody Information

Catalog Number: 86064-1-PBS	Clone ID: 242800G10	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human, mouse	Full name: within bgcn homolog (Drosophila)
Isotype: IgG	GenBank: BC006135	Gene ID: 84305
Purification Method: Protein A purification	Immunogen Catalog Number: Ag6276	

Detection Antibody Information

Catalog Number: 86064-2-PBS	Clone ID: 242800C3	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: within bgcn homolog (Drosophila)
Isotype: IgG	GenBank: BC006135	Gene ID: 84305
Purification Method: Protein A purification	Immunogen Catalog Number: Ag6276	

Applications

Tested Applications: Cytometric bead array	Range: 1.562-200 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

Product Information

MP02733-1 targets WIBG in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: WIBG Recombinant monoclonal antibody, PBS Only (Capture) 86064-1-PBS (242800G10). 100 µg. Concentration 1 mg/mL.

Detection antibody: WIBG Recombinant monoclonal antibody, PBS Only (Detector) 86064-2-PBS (242800C3). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP02733-1, WIBG Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86064-1-PBS. Detection antibody: 86064-2-PBS. Standard: Ag6276. Range: 1.562-200 ng/mL.



The mean WIBG concentration was determined to be 112.9 ng/mL in U2OS cell extract based on a 1.2 mg/mL extract load, 201.7 ng/mL in A549 cell extract based on a 1.2 mg/mL extract load.