

Mouse FAM3D Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP02761-1**

Capture Antibody Information

Catalog Number:
86090-4-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
250561F6
Reactivity:
mouse
GenBank:
NM_146050.2
Immunogen Catalog Number:
Eg2694

Conjugate:
Unconjugated
Full name:
oncoprotein induced transcript 1
Gene ID:
18300

Detection Antibody Information

Catalog Number:
86090-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
250561F4
Reactivity:
mouse, rat
GenBank:
NM_146050.2
Immunogen Catalog Number:
Eg2694

Conjugate:
Unconjugated
Full name:
oncoprotein induced transcript 1
Gene ID:
18300

Applications

Tested Applications:
Cytometric bead array

Range:
0.312-40 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP02761-1 targets FAM3D in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Mouse FAM3D Recombinant monoclonal antibody, PBS Only (Capture) 86090-4-PBS (250561F6). 100 µg. Concentration 1 mg/mL.

Detection antibody: Mouse FAM3D Recombinant monoclonal antibody, PBS Only (Detector) 86090-1-PBS (250561F4). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP02761-1, MOUSE FAM3D Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86090-4-PBS. Detection antibody: 86090-1-PBS. Standard: Eg2694. Range: 0.312-40 ng/mL



The mean FAM3D concentration was determined to be 2291.8 ng/mL in Mouse large intestine tissue extract based on a 4.5 mg/mL extract load.