

For Research Use Only



# TAC3 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP02770-1

## Capture Antibody Information

**Catalog Number:**  
86914-2-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
250567H3  
**Reactivity:**  
human  
**GenBank:**  
NM\_013251.3  
**Immunogen Catalog Number:**  
Eg2781

**Conjugate:**  
Unconjugated  
**Full name:**  
tachykinin 3  
**Gene ID:**  
6866

## Detection Antibody Information

**Catalog Number:**  
86914-1-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
250567E11  
**Reactivity:**  
human  
**GenBank:**  
NM\_013251.3  
**Immunogen Catalog Number:**  
Eg2781

**Conjugate:**  
Unconjugated  
**Full name:**  
tachykinin 3  
**Gene ID:**  
6866

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.312-40 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP02770-1 targets TAC3 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TAC3 Recombinant monoclonal antibody, PBS Only (Capture) 86914-2-PBS (250567H3). 100  $\mu$ g. Concentration 1 mg/mL.

Detection antibody: TAC3 Recombinant monoclonal antibody, PBS Only (Detector) 86914-1-PBS (250567E11). 100  $\mu$ g. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
**Storage buffer:**  
PBS only

For technical support and original validation data for this product please contact:



T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP02770-1, TAC3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86914-2-PBS. Detection antibody: 86914-1-PBS. Standard: Eg2781. Range: 0.312-40 ng/mL

The mean TAC3 concentration was determined to be 4.8 ng/mL in Human placenta tissue extract based on a 3.2 mg/mL extract load.