

# COX5A Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP02785-1**

## Capture Antibody Information

**Catalog Number:**  
86938-3-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
251697A7  
**Reactivity:**  
human  
**GenBank:**  
BC024240  
**Immunogen Catalog Number:**  
Ag2005

**Conjugate:**  
Unconjugated  
**Full name:**  
cytochrome c oxidase subunit Va  
**Gene ID:**  
9377

## Detection Antibody Information

**Catalog Number:**  
86938-2-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
251697H6  
**Reactivity:**  
human  
**GenBank:**  
BC024240  
**Immunogen Catalog Number:**  
Ag2005

**Conjugate:**  
Unconjugated  
**Full name:**  
cytochrome c oxidase subunit Va  
**Gene ID:**  
9377

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.312-40 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP02785-1 targets COX5A in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: COX5A Recombinant monoclonal antibody, PBS Only (Capture) 86938-3-PBS (251697A7). 100 µg. Concentration 1 mg/mL.

Detection antibody: COX5A Recombinant monoclonal antibody, PBS Only (Detector) 86938-2-PBS (251697H6). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
**Storage buffer:**  
PBS only

## Selected Validation Data



Cytometric bead array standard curve of MP02785-1, COX5A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86938-3-PBS. Detection antibody: 86938-2-PBS. Standard: Ag2005. Range: 0.312-40 ng/mL



The mean COX5A concentration was determined to be 505.7 ng/mL in mouse brain tissue extract based on a 6.3 mg/mL extract load.