

BRD2 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP02788-1**

Capture Antibody Information

Catalog Number: 85990-4-PBS	Clone ID: 250453G1	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: bromodomain containing 2
Isotype: IgG	GenBank: BC063840	Gene ID: 6046
Purification Method: Protein A purification	Immunogen Catalog Number: Ag17589	

Detection Antibody Information

Catalog Number: 85990-3-PBS	Clone ID: 250453A12	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: bromodomain containing 2
Isotype: IgG	GenBank: BC063840	Gene ID: 6046
Purification Method: Protein A purification	Immunogen Catalog Number: Ag17589	

Applications

Tested Applications: Cytometric bead array	Range: 1.25-40 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

Product Information

MP02788-1 targets BRD2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: BRD2 Recombinant monoclonal antibody, PBS Only (Capture) 85990-4-PBS (250453G1). 100 µg. Concentration 1 mg/mL.

Detection antibody: BRD2 Recombinant monoclonal antibody, PBS Only (Detector) 85990-3-PBS (250453A12). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP02788-1, BRD2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85990-4-PBS. Detection antibody: 85990-3-PBS. Standard: Ag17589. Range: 1.25-40 ng/mL.



The mean BRD2 concentration was determined to be 23.5 ng/mL in HeLa cell extract based on a 1.2 mg/mL extract load, 47.6 ng/mL in Jurkat cell extract based on a 1.2 mg/mL extract load, 21.7 ng/mL in A549 cell extract based on a 1.2 mg/mL extract load.