

MNK1 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP02797-1**

Capture Antibody Information

Catalog Number: 86966-3-PBS	Clone ID: 252096C8	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: MAP kinase interacting serine/threonine kinase 1
Isotype: IgG	GenBank: BC002755	Gene ID: 8569
Purification Method: Protein A purification	Immunogen Catalog Number: Ag0184	

Detection Antibody Information

Catalog Number: 86966-2-PBS	Clone ID: 252096B12	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: MAP kinase interacting serine/threonine kinase 1
Isotype: IgG	GenBank: BC002755	Gene ID: 8569
Purification Method: Protein A purification	Immunogen Catalog Number: Ag0184	

Applications

Tested Applications: Cytometric bead array	Range: 0.312-40 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	---	---

Product Information

MP02797-1 targets MNK1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: MNK1 Recombinant monoclonal antibody, PBS Only (Capture) 86966-3-PBS (252096C8). 100 µg. Concentration 1 mg/mL.

Detection antibody: MNK1 Recombinant monoclonal antibody, PBS Only (Detector) 86966-2-PBS (252096B12). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP02797-1, MNK1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86966-3-PBS. Detection antibody: 86966-2-PBS. Standard: Ag0184. Range: 0.312-40 ng/mL



The mean MNK1 concentration was determined to be 8.4 ng/mL in HeLa cell extract based on a 1.2 mg/mL extract load, 9.4 ng/mL in HEK-293 cell extract based on a 1.2 mg/mL extract load, 2.5 ng/mL in HepG2 cell extract based on a 1.2 mg/mL extract load.