

For Research Use Only

# PCNP Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP02799-1

## Capture Antibody Information

**Catalog Number:**  
86075-3-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
250576E5  
**Reactivity:**  
human  
**GenBank:**  
BC013916  
**Immunogen Catalog Number:**  
Ag30951

**Conjugate:**  
Unconjugated  
**Full name:**  
PEST proteolytic signal containing nuclear protein  
**Gene ID:**  
57092

## Detection Antibody Information

**Catalog Number:**  
86075-2-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
250576E7  
**Reactivity:**  
human  
**GenBank:**  
BC013916  
**Immunogen Catalog Number:**  
Ag30951

**Conjugate:**  
Unconjugated  
**Full name:**  
PEST proteolytic signal containing nuclear protein  
**Gene ID:**  
57092

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.938-60 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP02799-1 targets PCNP in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: PCNP Recombinant monoclonal antibody, PBS Only (Capture) 86075-3-PBS (250576E5). 100  $\mu$ g. Concentration 1 mg/mL.

Detection antibody: PCNP Recombinant monoclonal antibody, PBS Only (Detector) 86075-2-PBS (250576E7). 100  $\mu$ g. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
**Storage buffer:**  
PBS only

For technical support and original validation data for this product please contact:



T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MPO2799-1, PCNP Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86075-3-PBS. Detection antibody: 86075-2-PBS. Standard: Ag30951. Range: 0.938-60 ng/mL

The mean PCNP concentration was determined to be 25.9 ng/mL in A2780 cell extract based on a 1.2 mg/mL extract load, 38.9 ng/mL in Jurkat cell extract based on a 1.2 mg/mL extract load.