

# NUSAP1 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP02838-2

## Capture Antibody Information

**Catalog Number:**  
87034-5-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
252131D7  
**Reactivity:**  
human  
**GenBank:**  
BC001308  
**Immunogen Catalog Number:**  
Ag2654

**Conjugate:**  
Unconjugated  
**Full name:**  
nucleolar and spindle associated protein 1  
**Gene ID:**  
51203

## Detection Antibody Information

**Catalog Number:**  
87034-4-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
252131H3  
**Reactivity:**  
human  
**GenBank:**  
BC001308  
**Immunogen Catalog Number:**  
Ag2654

**Conjugate:**  
Unconjugated  
**Full name:**  
nucleolar and spindle associated protein 1  
**Gene ID:**  
51203

## Applications

**Tested Applications:**  
Sandwich ELISA

**Range:**  
46.9-1500 pg/mL (Sandwich ELISA)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP02838-2 targets NUSAP1 in immunoassays as a matched antibody pair. Validated in Sandwich ELISA.

Capture antibody: NUSAP1 Recombinant monoclonal antibody, PBS Only (Capture) 87034-5-PBS (252131D7). 100 µg. Concentration 1 mg/mL.

Detection antibody: NUSAP1 Recombinant monoclonal antibody, PBS Only (Detector) 87034-4-PBS (252131H3). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
**Storage buffer:**  
PBS only

## Selected Validation Data



Sandwich ELISA standard curve of MP02838-2, Human NUSAP1 Recombinant Matched Antibody Pair - PBS only. 87034-5-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag2654. 87034-4-PBS was HRP conjugated as the detection antibody. Range: 46.9-1500 pg/mL



The mean NUSAP1 concentration was determined to be 9.4 ng/mL in HEK-293 cell extract based on a 1.2 mg/mL extract load, 19.9 ng/mL in Jurkat cell extract based on a 1.2 mg/mL extract load and 0.2 ng/mL in A549 cell extract based on a 1.2 mg/mL extract load.