

For Research Use Only

Nidogen 2 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP02894-1

Capture Antibody Information

Catalog Number:
85299-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
252350B8
Reactivity:
human
GenBank:
BC035608
Immunogen Catalog Number:
Ag4643

Conjugate:
Unconjugated
Full name:
nidogen 2 (osteonidogen)
Gene ID:
22795

Detection Antibody Information

Catalog Number:
85299-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
252350B10
Reactivity:
human
GenBank:
BC035608
Immunogen Catalog Number:
Ag4643

Conjugate:
Unconjugated
Full name:
nidogen 2 (osteonidogen)
Gene ID:
22795

Applications

Tested Applications:
Cytometric bead array

Range:
1.25-80 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP02894-1 targets Nidogen 2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: NID2 Recombinant monoclonal antibody, PBS Only (Capture) 85299-3-PBS (252350B8). 100 μ g. Concentration 1 mg/mL.

Detection antibody: NID2 Recombinant monoclonal antibody, PBS Only (Detector) 85299-2-PBS (252350B10). 100 μ g. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:



T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MPO2894-1, Nidogen 2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85299-3-PBS. Detection antibody: 85299-2-PBS. Standard: Ag4643. Range: 1.25-80 ng/mL.

The mean Nidogen 2 concentration was determined to be 49.3 ng/mL in human placenta cell extract based on a 3.2 mg/mL extract load, HT-29 cell extract were served as a negative control.